

# Need something custom made? We Customize!

We can help transform your concept into reality using our closed cell foam and vinyl coating. Whether you are looking for a few items or thousands, email us with your design for a quote at [sales@danmarproducts.com](mailto:sales@danmarproducts.com).



## Hand Covers (6785)

Our all-foam Hand Covers are specially designed to help absorb and deflect the impact when an individual strikes out. They can be applied quickly with Velcro straps when aggressive and/or self-abusive behavior escalates, protecting both the wearer from injuring themselves on hard surfaces as well as possibly causing injury to person(s) nearby. Ventilated openings add comfort, and our soft, yet shock absorbent foam is easily cleaned with warm soapy water. Sold as a pair.



## Hard Shell Helmet (9821)

You can now provide individuals who are intent on banging their heads with the proper protection they need. These specially made protective helmets are comfortable to wear and effective in reducing the risk of head injuries. Its air ventilation holes allow for full ventilation without compromising head coverage. A hard polyethylene shell covers our shock absorbent foam lining that retains compression resistance and will not break down with multiple impacts. The adjustable chinstrap has a secure snap closure making it easy on and off. Many modifications are available to suit your needs, please call or email us for more information.



## Hard shell helmet with Spit mask

Stay Protected when confronted with a hostile or intoxicated member of the public or disturbed hospital patient or inmate; the risk of being spat on is real and protection from such a threat is essential. The hard shell with foam lining protects the individual from self-harm.

221 Jackson Industrial Drive Ann Arbor, MI 48103  
(800) 783-1998 (734) 761-8977 fax  
[www.danmarproducts.com](http://www.danmarproducts.com)



**DANMAR DEFENSE**  
PROTECTIVE EQUIPMENT



**800-783-1998**

[sales@danmarproducts.com](mailto:sales@danmarproducts.com)  
[www.danmarproducts.com](http://www.danmarproducts.com)



### Hard Shell Helmet with Face Guard and padded shell.

Our hard shell helmet has been padded on the outside to help protect others from injury from head butts. The protective face guard can be modified to protect against biting or spitting, and our custom strapping system with metal D-rings and clips makes it difficult to remove.




### Hard Shell Helmet with Face Guard (9822)

When you need full coverage, this model offers the best features. A clear, polycarbonate face guard is attached to a lightweight helmet, which is made from shock-absorbent foam inside a hard, ventilated, plastic shell. It has an eye level opening (1" X 6") for unobstructed vision and a mouth opening (1" x 3 3/4") for maximum ventilation. We also offer other variations for the faceguard to maximize coverage, for self-protection and for those around them. The bottom edge is lined with a strip of soft foam to help protect the chin. The foam is coated in vinyl, so it's easy to clean using just soap and water. We feature a wide range of custom options which we are happy to discuss with you.



### Hard Shell Helmet with Face Guard modifications.

Our Hard Shell Helmet with Face Guard has been modified to protect the user and those around them. The face guard features small ventilation holes that a pencil will not pierce which helps with self-abusiveness as well as protection from body fluids being expelled. Hard ear covers protect the ears and a rear hard shell extension and custom strapping system aids in preventing helmet removal.



### Helmet Modifications

3463 PERMANENT STABILIZING SIDE TABS

3459 HARD SHELL EAR COVERINGS



3470 EYEGLASS RELIEF AREA



3472 CUSTOM STRAPPING SYSTEM



3465 REAR HARD SHELL EXTENSION



### Hard Shell Helmet with Hockey Face Guard

Adding the hockey face guard to our hard helmet is the most durable option for facial protection. The added stabilizing side tabs add extra stabilization while also allowing the face guard being easily pivoted up and out of the face when needed. The custom strapping system is added to help prevent helmet removal.



### Hard Shell Comfy Cap (9709)

This 11.5 oz. helmet is our lightest Hard Shell made of a compression-molded inner foam liner attached to a rigid plastic exterior. The Hard Shell Comfy Cap is extremely well-ventilated and retains limited body heat.



### Comfy Cap (9708)

Our lightest, coolest helmet yet! This soft, lightweight helmet is made from compression-molded foam and weighs just 3 ounces. The Soft Comfy Cap provides maximum ventilation and retains limited body heat.



**DANMAR DEFENSE**  
PROTECTIVE EQUIPMENT

**800-783-1998**